

INFANT RADIANT WARMER NWS 101

WARM • FLEXIBLE • FRIENDLY



PHOENIX
Phoenix Medical Systems (P) Ltd.

CONTACT US NOW!

Authorised Representative and Distributor By

011 - 2417 3200 | Phone : +603 3324 8435

www.remedicss.com | remedic8@gmail.com

NO. 32-G, JALAN BAYU TINGGI 2A/KS6, TAMAN BAYU TINGGI, 41200 KLANG, SELANGOR



RADIATING warmth FRIENDLY to handle

The **NWS 101** provides radiant warmth for infants suffering from severe heat loss. It is suitable for all environments from labour rooms to neonatal ICUs.

FEATURES

Infrared (IR) radiation is beamed from a quartz heater by a parabolic reflector. This provides uniform warmth rapidly all over the bed. The heater can be operated in servo and manual modes.

The heater unit can be swivelled to accommodate X-ray units.

All the castors are provided with foot-operated brakes. The **NWS 101** may be shifted easily by just one person using the convenient handles.

The infant temperature is measured quickly using a high-sensitivity thermistor probe and displayed on a large, easy-to-read LED display. The display has a battery back-up and indicates the temperature even during power failures.

The control panel is ergonomically designed with soft-touch keys. This modular panel may be easily removed if servicing is required. An easily accessible quick reference is provided on the side of the **NWS 101**

The **NWS 101** is equipped with a dazzle-free observation light for examination of the infant.

Height adjustable IV pole is provided. Accessories like Infusion pump can be attached on the IV Pole

Storage space provided to place the baby things

The **NWS 101** is backed by the many years of experience **PHOENIX** has in designing and supplying sophisticated neonatal care equipment worldwide.



PHOENIX endeavours to meet your every need. If your practice has any special requirement, call **PHOENIX**. We will be happy to provide you a solution.

CE 0123

Call No: 7.2.3 - 28 / 04 / June 2020



Authorised Representative and Distributor By
REMEDIC SERVICES & SUPPLIES SDN BHD
NO. 32-G, JALAN BAYU TINGGI 2A/KS6, TAMAN BAYU TINGGI, 41200 KLANG
SELANGOR, MALAYSIA

☎ 011 - 2417 3200 | Phone : +603 3324 8435
✉ www.remedicss.com | remedic8@gmail.com